The SpecialNeeds™ Feeder

The Mini- SpecialNeeds™ Feeder

Manufactured by Medela, Inc.
P.O. Box 660, McHenry, IL 60051-0660
1101 Corporate Dr., McHenry, IL 60050
Phone: 800-435-8316 or 815-363-1166 Fax: 815-363-1246
Email: customer.service@medela.com

Medela Canada, Inc.
4090B Sladeview Crescent, Unit 2
Mississauga, Ontario L5L 5Y5
Phone: 905-608-7272 or 800-435-8316 Fax: 905-608-8720
Email: info@medela.ca

www.medela.com

Medela is a registered trademark and SpecialNeeds and
Mini-SpecialNeeds are trademarks of Medela
Need Translation.
Need Translation.

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
**SpecialNeeds™
Feeder**

Table of Contents

How Does the SpecialNeeds™ Feeder Help	1
Cleaning	4
Using the SpecialNeeds or Mini-SpecialNeeds Feeder	5
Feeding Your Baby	8
Helpful Hints	9
Replacement of parts	12

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Dear Parent,

Does your baby have difficulty feeding? Be reassured, you are not alone. Many babies born each year experience feeding problems for a wide variety of reasons. The SpecialNeeds™ Feeder has been developed specifically to help your baby. The Mini-SpecialNeeds™ Feeder is a smaller version for smaller or premature infants.

The SpecialNeeds feeder was invented by a mother whose baby was born with Stickler's Syndrome. She experienced the misery of feeding problems and recognized the need for a specialized feeder. Five years of research and development have resulted in a product which really can transform feeding from a traumatic experience into a relaxed and rewarding activity for both you and your baby. Its innovative design provides the answer to lengthy, distressful feedings and an alternative to nasogastric tube or spoon feedings. It is suitable for all babies for whom oral feeding is recommended.

First introduced to the United Kingdom in 1987, the SpecialNeeds Feeder has been widely used for babies suffering from ailments such as cleft palate, chromosomal abnormality, cardiac conditions, neuromuscular weakness, mental handicap, cerebral palsy and other conditions.

**Whatever your baby's condition, the
SpecialNeeds Feeder is here to help you.**

How Does the SpecialNeeds™ Feeder Help?

Rewards the weakest suck

In conventional feeders, much of the baby's effort is wasted on compression of air within the container and movement of feed within the system. However, the unique SpecialNeeds teat is separated from the container by a one-way valve. Once filled, milk is held in the teat and cannot flow back into the container. All your baby's sucking effort is directed onto the contents of the teat, so even the weakest suck gets results.

Never floods

Enlarging the hole of an ordinary teat does not solve your feeding problems. Your baby cannot stop the flow of feed in order to relax and breathe. Panic and choking soon result. The SpecialNeeds Feeder has a slit-valve in the mouthpiece which only opens when your baby sucks. Your baby has complete control and is never overwhelmed. It also has the advantage of variable flow. Graduation from a large opening, down to a shut position, can be achieved to suit the individual needs of your baby. Whichever opening you select, it will always shut between sucks.

May reduce gas, vomiting and colic


Many babies swallow air while using some ordinary feeders. This is because a vacuum builds up inside the system and the baby must release his lips from the teat to allow air back in. Consequently, air is swallowed with the next mouthful of feed causing gas, vomiting and colic.

The SpecialNeeds™ Feeder may reduce air-swallowing. It has a special air inlet groove which allows air in at the same rate as feed is released from the system. There is no vacuum to work against your baby's sucking efforts and air-swallowing is avoided. Even babies with cleft palates benefit significantly from this feature.

You can help your baby to feed

If your baby needs some additional help, you can gently squeeze and release the reservoir of the teat. This will deliver a controlled amount of feed into your baby's mouth. The full teat has a natural resistance which prevents you from squeezing too hard.

The SpecialNeeds™ and Mini-SpecialNeeds™ Feeder

Cleaning

After each feeding: Flush cold water through the slit in the teat to remove as much feed as possible. (Do not turn the teat inside out or push brushes, etc. through the slit.) Take the feeder apart, making sure that the valve membrane is separated from the disc. Wash all the parts in warm soapy water. Rinse thoroughly with cold water.

Each day: All disassembled parts of the feeder should be sanitized by boiling, preferably in distilled water, for 20 minutes. If hard water deposits occur use distilled water or occasionally boil in a mixture of 3 parts water and 1 part vinegar. Parts can be autoclaved to a maximum of 272°F for 3 minutes at 29 PSIG or 250°F for 15 minutes at 15 PSIG. To avoid damage during sterilization, sterilize parts without applying external pressure on parts.

Using The SpecialNeeds™ or Mini-SpecialNeeds™ Feeder

IMPORTANT

Consult your doctor before using the SpecialNeeds or Mini-SpecialNeeds Feeder for the first time, in order to make sure that oral feeding is recommended for your baby. Your doctor should tell you the appropriate amount of weight your baby should gain in a specified time frame. If your baby does not gain this amount of weight, contact your doctor or other health care provider. Always inspect parts thoroughly before each use.


Components

Your SpecialNeeds or Mini-SpecialNeeds Feeder should include:

- 1-Container, 1-Teat, 1-Disc,
- 1- Valve membrane, 1-Collar.

Assembly

1. Press valve membrane onto the upper side of the disc, so that the stud goes completely through the center hole.

2. Fill the container with room temperature breastmilk or supplement recommended by physician to the level that the baby is presently taking.

3. Slip the teat into the collar.


4. Put assembled valve membrane into the teat. Make sure that the valve membrane and high rim of the disc are facing the inside of the teat.

5. Place the teat assembly over the container and screw all the parts together using the collar.

Filling the teat

The feeder is designed for use with room temperature feed and should not be used with hot liquids.

1. Hold the feeder upright and **SQUEEZE** some air out of the teat.
2. **KEEP SQUEEZING** and tip the feeder upside down.
3. Then **RELEASE** and some feed will enter the teat.
4. **REPEAT** steps 1-3 several times until the teat is almost full. As your baby feeds, the teat will automatically refill.



Controlling the flow

The mouthpiece of the teat has an opening with a slit-valve, which can be adjusted to suit your baby. The flow rate is shown by the three lines on the teat—longest line for maximum flow, medium line for medium flow, shortest line for zero flow.

How the slit valve works

- a) Slit horizontal: pressure keeps valve closed for zero flow
- b) Slit midway: valve partly opens for medium flow
- c) Slit vertical: pressure opens valve for maximum flow


Measuring your baby's intake

The whole teat holds approx.


The mouthpiece holds approx.

	Mini-SpecialNeeds	SpecialNeeds
The whole teat holds approx.	27ml	25ml
The mouthpiece holds approx.	5ml	4ml


The contents of the container can be read on the scale on the side of the container.

Standard SpecialNeeds Teat

Mini-SpecialNeeds Teat

Feeding Your Baby

1. Relax

2. Select the flow rate

To adjust the rate of flow, turn the feeder so that the required line (zero, medium or maximum) points toward your baby's nose. If baby has a weak suck, start with the teat positioned between the medium and longest line, then adjust as necessary. Use the zero flow if your baby needs to get used to the teat before feeding begins. When your baby has settled, gently increase the flow. Reduce it and start again if your baby shows signs of distress.

3. Helping your baby to feed

If your baby needs help, the following technique may be useful. However, before using this feature for the first time, we recommend that you practice squeezing some water into the palm of your hand. Hold the feeder with the collar resting in the angle between your thumb and first finger. Gently squeeze and release the straight body section, or reservoir, of the teat. This technique will put a small amount of fluid into your baby's mouth. Always give your baby time to swallow before squeezing again. Try to establish a good rhythm in harmony with your baby's suck and swallow. You will soon feel confident in your ability to feed your baby.

NOTE: *The amount given by squeezing is not affected by the variable flow positions of the teat.*

Helpful Hints

Teat not refilling properly

If you are using thickened feed or additives and the teat becomes floppy, unscrew the collar and check that the air-groove in the disc has not become blocked by sediment. You may need to use finger pressure to help the flow of feed. The teat will also not fill properly if the disc is upside down.

Bubbles of air in teat

If air is entering the teat as you feed, tilt the container a little more so that feed completely covers the valve membrane and disc. Sometimes on maximum flow, the baby's sucking action may cause an occasional bubble of air to enter via the mouthpiece. If this happens often, set the flow to just below maximum.

When to use the SpecialNeeds™ or Mini-SpecialNeeds™

The Mini-SpecialNeeds was designed for smaller or premature infants. Consult with your physician or other health care provider for advice on which product to use for your baby.

Feeder leaking

If the feeder leaks, check that the feeder has been correctly assembled. The stud of the valve membrane must be properly inserted through the hole in the disc. The collar must be firmly screwed onto the container. Make sure the feed is not too hot. The feeder is designed for use with room temperature feed and should not be used with hot liquids.

Check to see if either the valve membrane or teat needs to be replaced. Also, if using a container other than the Medela 80ml or 150ml containers, some leakage may occur around the collar and container connection.

Feeder not working

Check that the slit in the mouthpiece is opening properly. Hold it vertically between your fingers and squeeze it open. The slit should be about 3mm long. Sometimes, after storage, the slit may close up. If this occurs, apply warm soapy water and gently ease it open. Thoroughly rinse in cold water before use.

Traveling

To keep the cost of the unit as low as possible, we do not provide a travel lid. Instead, invert the teat and disc so that the teat is inside the container, then screw on the collar. Alternatively, put the feeder inside a sanitized, water-tight plastic container.

Using expressed breast milk

For mothers feeding expressed breast milk, please note that the container of the SpecialNeeds Feeder can be directly attached to the breastshield of any Medela Breastpump Accessory Kit.

Container size

The SpecialNeeds Feeder is available with either an 80ml or 150ml container. The Mini-SpecialNeeds Feeder is available with an 80ml container. Other containers may be used, however you may experience some leakage around the collar and the container.

Baby rejects teat

If your baby rejects the teat, squeeze a little feed onto your finger and let him taste it.

Touching your baby's cheek with the teat may encourage him to feed by stimulating the natural rooting reflex.

Remember to be relaxed and patient. The initial success of using this device may not be measured in a large volume of intake, but rather by your baby learning to coordinate a routine of sucking, swallowing and breathing.

Baby still not feeding well

Ask your Speech Therapist, Lactation Consultant, Physician or other health care professional for advice. Feeding can often be improved by adjusting the baby's position and/or simple stimulation techniques.

Removing Teat

If there is some feed left in the feeder, hold it horizontally over a sink before unscrewing the collar. The disc should stay in the teat, preventing feed from spilling out. Quickly stand the container upright.

Replacement of Parts

As a guide, we suggest that the teat and valve membrane are replaced every five weeks.

Parts are worn out if:

- The slit in the teat has grown longer than 4mm.
- The valve membrane is warped and no longer sits flat on the disc.

Always inspect parts thoroughly before each use.

Parts may also be ordered by calling 1-800-435-8316 and ordering with Discover Card, Mastercard, Visa, or American Express.

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© 2006 Medela, Inc.