



neonatal solutions

# Dual Cap

U.S. Design Patent D812,222

The Medela Dual Cap provides a convenient way for your pharmacy to continue the use of traditional oral slip-tip syringes while allowing the NICU to benefit from the transition to the ENFit® enteral connectors. Medela's Dual Cap does not require twisting to apply the cap to the syringe. This feature results in one less step during filling, which may result in saved time and labor.

## Benefits:


- Pharmacy can use the same syringes and filling methods
- Consists of a 2-piece cap for oral slip-tip syringes
- Removal of both components of the Dual Cap allows PO dosing identical to current use
- Removal of first cap facilitates connection to ENFit feeding tubes
- Designed with Low Dose technology, which provides accurate dosing through the feeding tube

## Product Specifications:


Item #	Description	Sterile or Non-Sterile	Qty/Case (each)
ENFDUALCAPLD	ENFit to TwistLok® connector, female to female with cap	Sterile	500 (50 packs of 10) per case
ENFDUALCAPNSLD	ENFit to TwistLok connector, female to female with cap	Non-Sterile	500 caps per pouch
ENFDUALCAPCLD	ENFit LDT to TwistLok connector, female to female with cap for Covidien oral syringe.	Sterile	500 (50 packs of 10) per case
ENFDUALCAPCNSLD	ENFit LDT to TwistLok connector, female to female with cap for Covidien oral syringe.	Non-Sterile	500 caps per pouch

To order or for more information, contact your Medela Sales Consultant.


## Just 3 easy steps....



**1** No change to the PO delivery process.



**2** Removeable cap for connection to the ENFit feeding tube.



**3** Low-dose connector provides accurate dosing.

Content appearance may vary from picture.

Medela wordmark and logo and TwistLok are registered in the U.S. Patent and Trademark Office. ENFit is a federally registered trademark of GEDSA and is used with permission. Neonatal Solutions logo is a trademark of Medela.  
1549104 A 318 © 2018 Medela 4550174