

BREASTPUMP TUBING

Reference Sheet

Updated 03/09



Item



#8007213

(Replaces item #8007179 and #8007194)

42" PVC Tubing that features one two-tone yellow adapter and one pale yellow adapter.

For use with the Classic™, Lactina®, Symphony®
Not recommend to be used with Pump In Style Advanced.



#87214

(Replaces item #8007180 & #8007214)

42" Silicone Tubing set. Each tube features a two-tone yellow adapter and one pale yellow adapter.

For use with the Classic™, Lactina®, Symphony® and any Pump in Style® Advanced sold prior to August 2006 (old faceplate)



#87212

(Replaces item #8007156 & #8007212)

42" PVC Tubing set. Each tube features pale yellow adapter and one soft end.

For use with all PNSO** and any PNSA sold after July 2006 (new faceplate)

Includes all Limited Edition** & Metro Styles



*Not pictured



#8087018

24" Silicone tubing with one soft end and one white adapter

For use with item #67274- DoubleSelect™



#8007215

30" PVC tubing featuring two yellow adapters

For use with item #67050 - Swing™



** = Discontinued product